

# Anti-SCGN / Secretagogin antibody



Catalog Number: 100713

## Product name

Anti-SCGN / Secretagogin antibody

## Immunogen

[Human SCGN / Secretagogin \(His Tag\) recombinant protein](#)

## Specificity

Human SCGN / Secretagogin

## Antibody description

Rabbit monoclonal to SCGN / Secretagogin

## Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human SCGN / Secretagogin (rh SCGN / Secretagogin; NP\_008929.2; Asp2-Pro276).

## Formulation

0.2 µm filtered solution in PBS with 5% trehalose

## Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid

contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

IHC-P

## Dilutions

IHC-P: 1-10 µg/mL

## Validations



SCGN / Secretagogin Antibody, Rabbit MAb, Immunohistochemistry

Immunochemical staining of human SCGN in human pancreas (from 2 donors) with rabbit monoclonal antibody (5 µg/mL, formalin-fixed paraffin embedded sections).