

# FOXQ1 antibody



Catalog Number: 110713

## Product name

FOXQ1 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

FOXQ1 Rabbit Polyclonal antibody. Positive WB detected in mouse kidney tissue, rat kidney tissue. Observed molecular weight by Western-blot: 42 kDa

## Preparation

This antibody was obtained by immunization of FOXQ1 recombinant protein (Accession Number: NM\_033260). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



mouse kidney tissue were subjected to SDS PAGE followed by western blot with Catalog No:110713(FOXQ1 antibody) at dilution of 1:800