

# Anti-TRAPPC2 Rabbit antibody



Catalog Number: 167713

## Product name

Anti-TRAPPC2 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

TRAPPC2 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human TRAPPC2

## Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,  
40%Glycerol, 0.01% sodium azide and 0.05% BSA

## Storage

Store at 4°C short term. Store at -20°C long term.  
Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of TRAPPC2 in HeLa,HL-60  
using TRAPPC2 Rabbit mAb(1:1000 diluted)