

Human UBE2G2 (GST tag) recombinant protein



Catalog Number: 519713

General Information

Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 G2, Ubiquitin Carrier Protein G2, Ubiquitin-Protein Ligase G2, UBE2G2

Protein Construction

A DNA sequence encoding the Human UBE2G2 Met1-Leu165 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human UBE2G2 (GST tag) recombinant protein