

EYA4 antibody

Catalog Number: 110423

Product name

EYA4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

EYA4 Rabbit Polyclonal antibody. Positive WB detected in mouse skeletal muscle tissue, HeLa cells, mouse heart tissue. Positive IF detected in HeLa cells. Observed molecular weight by Western-blot: 67 kDa

Preparation

This antibody was obtained by immunization of EYA4 recombinant protein (Accession Number: NM_172103). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:110423(EYA4 Antibody) at dilution of 1:1000



Immunofluorescent analysis of HeLa cells using Catalog No:110423(EYA4 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)