

# MAP4K1 antibody



Catalog Number: 112423

## Product name

MAP4K1 antibody

## Specificity

Human; other species not tested.

## Antibody description

MAP4K1 Rabbit Polyclonal antibody. Positive WB detected in Raji cells, Jurkat cells. Positive IP detected in Jurkat cells. Observed molecular weight by Western-blot: 90 kDa

## Preparation

This antibody was obtained by immunization of MAP4K1 recombinant protein (Accession Number: NM\_007181). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IP: 1:500-1:5000

## Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:112423(MAP4K1 antibody) at dilution of 1:1000

IP Result of anti-MAP4K1 (IP:Catalog No:112423, 4ug; Detection:Catalog No:112423 1:1000) with Jurkat cells lysate 3200ug.