

ORC6 antibody

Catalog Number: 113423

Product name

ORC6 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

ORC6 Rabbit Polyclonal antibody. Positive WB detected in Jurkat cells, MCF7 cells, mouse spleen tissue. Positive IHC detected in human cervical cancer tissue. Observed molecular weight by Western-blot: 28 kDa

Preparation

This antibody was obtained by immunization of ORC6 recombinant protein (Accession Number: NM_014321). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200

Validations



Jurkat cells were subjected to SDS PAGE followed by western blot with Catalog No:113423(ORC6L antibody) at dilution of 1:800



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:113423(ORC6L antibody) at dilution of 1:100 (under 10x lens)