

RAB24 antibody

Catalog Number: 114423

Product name

RAB24 antibody

Immunogen

[Human RAB24 Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse; other species not tested.

Antibody description

RAB24 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, HEK-293 cells, human skeletal muscle tissue. Positive IP detected in mouse brain tissue. Positive IHC detected in human gliomas tissue. Observed molecular weight by Western-blot: 23 kDa

Preparation

This antibody was obtained by immunization of RAB24 recombinant protein (Accession Number: NM_130781). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB, IP

Dilutions


Recommended Dilution:

WB: 1:200-1:1000


IP: 1:200-1:1000

IHC: 1:20-1:200

Validations




mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114423(RAB24 Antibody) at dilution of 1:300



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:114423(RAB24 antibody) at dilution of 1:50 (under 10x lens)

RAB24 antibody

Catalog Number: 114423



IP Result of anti-RAB24 (IP:Catalog No:114423, 3ug; Detection:Catalog No:114423 1:300) with mouse brain tissue lysate 4000ug.