

INHBC Polyclonal Antibody



Catalog Number: 165423

Product name

INHBC Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to INHBC

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 103-352 of human INHBC (NP_005529.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000


IHC 1:100 - 1:200

Validations



Western blot - INHBC Polyclonal Antibody


Western blot analysis of extracts of various cell lines, using INHBC antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure



time: 90s.

Immunohistochemistry - INHBC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat kidney using INHBC antibody at dilution of 1:100



(40x lens).


INHBC Polyclonal Antibody



Catalog Number: 165423

Immunohistochemistry - INHBC Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human colon using INHBC antibody at dilution of



1:100 (40x lens).

Immunohistochemistry - INHBC Polyclonal Antibody


Immunohistochemistry of paraffin-embedded human stomach using INHBC antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - INHBC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human esophagus using INHBC antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - INHBC Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using INHBC antibody at dilution of 1:100 (40x lens).