

# SMAD9 antibody



Catalog Number: 115423

## Product name

SMAD9 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

SMAD9 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293 cells, A549 cells, HepG2 cells, mouse lung tissue. Observed molecular weight by Western-blot: 52kd

## Preparation

This antibody was obtained by immunization of SMAD9 recombinant protein (Accession Number: NM\_001378621). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115423(SMAD9 antibody) at dilution of 1:500