

Phospho-LIMK1-T508 pAb

Catalog Number: 166423

Product name

Phospho-LIMK1-T508 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-LIMK1-T508

Preparation

Antigen: A phospho specific peptide corresponding to residues surrounding T508 of human LIMK1

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:100 - 1:200

Validations




Western blot - Phospho-LIMK1-T508 pAb

Western blot analysis of extracts from HeLa cells, using Phospho-LIMK1-T508 antibody .



Immunohistochemistry - Phospho-LIMK1-T508 pAb

Immunohistochemistry of paraffin-embedded human breast carcinoma using Phospho-LIMK1-T508 antibody .



Immunofluorescence - Phospho-LIMK1-T508 pAb

Immunofluorescence analysis of methanol-fixed HeLa cells using Phospho-LIMK1-T508 antibody .