

# MEST antibody



Catalog Number: 107423

## Product name

MEST antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

MEST Mouse Monoclonal antibody. Positive IHC detected in human prostate cancer tissue. Positive WB detected in PC-3 cells. Observed molecular weight by Western-blot: 34 kDa

## Preparation

This antibody was obtained by immunization of MEST recombinant protein (Accession Number: NM\_002402). Purification method: Protein A purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

IHC: 1:50-1:500

## Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107423(MEST antibody) at dilution of 1:1000