

Human ACYP2/ACYP (GST tag) recombinant protein



Catalog Number: 518423

General Information

Gene Name Synonym

Acylphosphatase-2, ACYP2, Acylphosphatase, muscle type isozyme, Acylphosphate phosphohydrolase 2, ACYP

Protein Construction

A DNA sequence encoding the Human ACYP2/ACYP Ser2-Tyr99 is expressed with a GST tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a 0.2 µm filtered solution of PBS, pH 7.4.

Usage Guide

Storage

Store at <-20°C, minimize freeze-thaw cycles.

SDS-PAGE



Human ACYP2/ACYP (GST tag) recombinant protein