

GADL1 antibody

Catalog Number: 110823

Product name

GADL1 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GADL1 Rabbit Polyclonal antibody. Positive IP detected in mouse brain tissue. Positive WB detected in mouse cerebellum tissue, human brain tissue. Observed molecular weight by Western-blot: 47 kDa

Preparation

This antibody was obtained by immunization of GADL1 recombinant protein (Accession Number: BC093701). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



mouse cerebellum tissue were subjected to SDS PAGE followed by western blot with Catalog No:110823(GADL1 antibody) at dilution of 1:300



IP Result of anti-GADL1 (IP:Catalog No:110823, 4ug; Detection:Catalog No:110823 1:300) with mouse brain tissue lysate 4000ug.