

Anti-ATG7 Rabbit antibody



Catalog Number: 166823

Product name

Anti-ATG7 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

ATG7 Rabbit polyclonal antibody

Preparation

Antigen: A synthetic peptide of human Apg7

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at -20°C. Stable for 12 months from date of receipt.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/2000-1/10000

ICC/IF: 1/20-1/100

Validations



Western blot detection of Apg7 in K562, Rat Brain, C6, 3T3, HeLa cell lysates using Apg7 Rabbit pAb(1:1000 diluted). Predicted band size: 78kDa. Observed band size: 78kDa.