

# ACAA1 antibody

Catalog Number: 107823

## Product name

ACAA1 antibody

WB: 1:500-1:5000

IHC: 1:20-1:200

## Specificity


Human, Mouse, Rat; other species not tested.

IF: 1:10-1:100

## Antibody description

ACAA1 Rabbit Polyclonal antibody. Positive WB detected in L02 cells, A549 cells, human lung tissue. Positive IF detected in HeLa cells. Positive IHC detected in human thyroid cancer tissue, human liver cancer tissue, human prostate cancer tissue. Observed molecular weight by Western-blot: 41 kDa

## Validations




## Preparation

This antibody was obtained by immunization of ACAA1 recombinant protein (Accession Number: NM\_001607). Purification method: Antigen affinity purified.

L02 cells were subjected to SDS PAGE followed by western blot with Catalog No:107823(ACAA1 antibody) at dilution of 1:1000

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

Immunohistochemistry of paraffin-embedded human thyroid cancer tissue slide using Catalog No:107823(ACAA1 Antibody) at dilution of 1:50 (under 10x lens)

## Ig Type


Rabbit IgG

## Applications

ELISA, WB, IF, IHC

## Dilutions

Recommended Dilution:




# ACAA1 antibody

Catalog Number: 107823

---

Immunohistochemistry of paraffin-embedded human thyroid cancer tissue slide using Catalog No:107823(ACAA1 Antibody) at dilution of 1:50 (under 40x lens)



Immunofluorescent analysis of HeLa cells using Catalog No:107823(ACAA1 Antibody) at dilution of 1:25 and Rhodamine-Goat anti-Rabbit IgG