

PPID antibody



Catalog Number: 114123

Product name

PPID antibody

Recommended Dilution:

WB: 1:500-1:5000

Specificity

Human, Mouse, Rat; other species not tested.


IP: 1:500-1:5000

IHC: 1:20-1:200

Antibody description

PPID Rabbit Polyclonal antibody. Positive FC detected in HeLa cells. Positive IHC detected in human cervical cancer tissue. Positive WB detected in mouse brain tissue, A549 cells, human brain tissue, Jurkat cells, MCF7 cells, NIH/3T3 cells, Raji cells, rat brain tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 38-41 kDa

Validations



Preparation

This antibody was obtained by immunization of PPID recombinant protein (Accession Number: NM_005038). Purification method: Antigen affinity purified.

mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:114123(PPID antibody) at dilution of 1:1000

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal


Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC, FC, IP


Dilutions



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:114123(PPID antibody) at dilution of 1:50 (under 10x lens)


PPID antibody

Catalog Number: 114123



Immunohistochemical of paraffin-embedded human cervical cancer using Catalog No:114123(PPID antibody) at dilution of 1:50 (under 40x lens)

IP Result of anti-PPID (IP:Catalog No:114123, 4ug; Detection:Catalog No:114123 1:1000) with mouse brain tissue lysate 4000ug.



1X10⁶ HeLa cells were stained with 0.2ug PPID antibody (Catalog No:114123, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L) with dilution 1:1500.

