

Anti-Cadherin 16 Mouse antibody



Catalog Number: 169123

Product name

Anti-Cadherin 16 Mouse antibody

Specificity

Human

Antibody description

Cadherin 16 Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant fragment of human CDH16 (AA: EX(558-721)) expressed in E. Coli.

Formulation

Purified antibody in PBS with 0.05% sodium azide.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

WB, IHC-P, ICC/IF, FC


Validations



Immunofluorescence analysis of A431 cells using CDH16 mouse mAb (green). Blue



Immunohistochemical analysis of paraffin-embedded kidney tissues using CDH16 mouse




mAb with DAB staining.

Western blot analysis using CDH16 mAb against HEK293 (1) and CDH16 (AA


Anti-Cadherin 16 Mouse antibody




Catalog Number: 169123



Flow cytometric analysis of A431 cells using CDH16 mouse mAb (green) and negative control (red).



Western blot analysis using CDH16 mAb against



human CDH16 (AA)

Western blot analysis using CDH16 mouse mAb against A431 (1), Jurkat (2) cell lysate.