

# HOP2 antibody

Catalog Number: 111523

## Product name

HOP2 antibody

IP: 1:200-1:2000

IHC: 1:20-1:200

## Specificity


Human; other species not tested.

IF: 1:10-1:100

## Antibody description

HOP2 Rabbit Polyclonal antibody. Positive WB detected in K-562 cells, Jurkat cells. Positive IP detected in Jurkat cells. Positive IHC detected in human lung cancer tissue. Positive IF detected in HepG2 cells. Observed molecular weight by Western-blot: 24 kDa and 29 kDa

## Validations




## Preparation

This antibody was obtained by immunization of HOP2 recombinant protein (Accession Number: NM\_013290). Purification method: Antigen affinity purified.

K-562 cells were subjected to SDS PAGE followed by western blot with Catalog No:111523(PSMC3IP antibody) at dilution of 1:500

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.



## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB, IHC, IF, IP

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

Immunohistochemical of paraffin-embedded human lung cancer using Catalog No:111523(PSMC3IP antibody) at dilution of 1:50 (under 10x lens)



# HOP2 antibody

Catalog Number: 111523

---

Immunofluorescent analysis of HepG2 cells, using PSMC3IP antibody Catalog No:111523 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit

IP Result of anti-PSMC3IP (IP:Catalog No:111523, 4ug; Detection:Catalog No:111523 1:500) with Jurkat cells lysate 3200ug.

