

CYP4A11 Polyclonal Antibody



Catalog Number: 164523

Product name

CYP4A11 Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to CYP4A11

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CYP4A11 (NP_000769.2).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100

Validations



Immunohistochemistry - CYP4A11 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse kidney using CYP4A11 antibody at dilution of 1:100 (40x lens).