

# Anti-Caspase-1 Rabbit antibody



Catalog Number: 167523

## Product name

Anti-Caspase-1 Rabbit antibody

## Specificity

Human, Rat

## Antibody description

Caspase-1 Rabbit polyclonal antibody

## Preparation

Antigen: A synthesized peptide derived from human Caspase-1

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC-P

## Dilutions


WB: 1/500-1/2000

IHC: 1/50-1/200

## Validations



Immunohistochemical analysis of paraffin-embedded human liver, using Caspase-1 Antibody.



Western blot analysis of Caspase-1 in HeLa cell lysate.