

SECTM1 antibody



Catalog Number: 107523

Product name

SECTM1 antibody

Specificity

Human; other species not tested.

Antibody description

SECTM1 Mouse Monoclonal antibody. Positive FC detected in K-562 cells. Positive WB detected in PC-3 cells, K-562 cells, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 27 kDa

Preparation

This antibody was obtained by immunization of SECTM1 recombinant protein (Accession Number: XM_005256392). Purification method: Caprylic acid/ammonium sulfate precipitation.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgM

Applications


ELISA, WB, FC

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:107523(SECTM1 Antibody) at dilution of 1:1000



1X10⁶ K-562 cells were stained with 0.2ug SECTM1 antibody (Catalog No:107523, red) and control antibody (blue). Fixed with 90% MeOH blocked with 3% BSA (30 min). Alexa Fluor 488-conjugated AffiniPure Goat Anti-Mouse IgG(H+L) with dilution 1:1000.