

# Anti-S100A10 (6F4) Mouse antibody



Catalog Number: 168523

## Product name

Anti-S100A10 (6F4) Mouse antibody

## Specificity

Human, Mouse, Rat

## Antibody description

S100A10 (6F4) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human S100A10 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a

## Applications

WB, IP, ICC/IF


## Dilutions

WB: 1/1000


;ICC: 1/200

;IP: 1/100


## Validations



Immunoprecipitation analysis of HeLa cell lysates using S100A10 mouse mAb.



Western blot detection of S100A10 in C6,3T3 and HeLa cell lysates using S100A10 mouse mAb (1:1000 diluted). Predicted band size:11KDa. Observed band size:11KDa.



Immunocytochemistry staining of HeLa cells fixed with 4% Paraformaldehyde and using anti-S100A10 mouse mAb (dilution 1:200).