

BAP1 Polyclonal Antibody



Catalog Number: 164923

Product name

BAP1 Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

IF 1:50 - 1:200

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to BAP1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 460-729 of human BAP1 (NP_004647.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC, IF


Dilutions

Validations



Immunohistochemistry - BAP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast using BAP1 antibody at dilution of



1:150 (40x lens).

Immunohistochemistry - BAP1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using BAP1 antibody at dilution of 1:150 (40x lens).