

# TDO2 antibody

Catalog Number: 115923

---

## Product name

TDO2 antibody

## Immunogen

[Human TDO2 Recombinant protein \(His tag\)](#)

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

TDO2 Rabbit Polyclonal antibody. Positive WB detected in mouse liver tissue, rat liver tissue. Observed molecular weight by Western-blot: 35 kDa, 48 kDa

## Preparation

This antibody was obtained by immunization of TDO2 recombinant protein (Accession Number: NM\_005651). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

ELISA, WB

## Dilutions

Recommended Dilution:  
WB: 1:200-1:2000

## Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115923(TDO2 Antibody) at dilution of 1:600