

Anti-TrkB Rabbit antibody



Catalog Number: 168923

Product name

Rabbit IgG

Anti-TrkB Rabbit antibody

Applications

Specificity

WB, IHC-P, ICC/IF

Human, Mouse, Rat

Dilutions

Antibody description

1:2000-1:10000

TrkB Rabbit monoclonal antibody

1:100

Preparation


1:50

Antigen: Recombinant protein of human TrkB

Validations

Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.



Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Western blot detection of TrkB (pan) in Rat Brain cell using TrkB (pan) Rabbit mAb(1:1000 diluted). Predicted band size:92kDa.Observed band size:120-140kDa.

Ig Type