

# DHRS2 antibody

Catalog Number: 109923

## Product name

DHRS2 antibody

## Specificity

Human; other species not tested.

## Antibody description

DHRS2 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in human liver tissue. Observed molecular weight by Western-blot: 30 kDa

## Preparation

This antibody was obtained by immunization of DHRS2 recombinant protein (Accession Number: BC007339). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:1000

## Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:109923(DHRS2 antibody) at dilution of 1:500



IP Result of anti-DHRS2 (IP:Catalog No:109923, 4ug; Detection:Catalog No:109923 1:300) with HeLa cells lysate 4000ug.