

Human NOL3/ARC/NOP (His tag) recombinant protein



Catalog Number: 519923

General Information

Gene Name Synonym

Nucleolar Protein 3, Apoptosis Repressor With CARD, Muscle-Enriched Cytoplasmic Protein, Myp, Nucleolar Protein of 30 kDa, Nop30

Protein Construction

A DNA sequence encoding the Human NOL3/ARC/NOP Met1-Ser208 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 25mM TrisHCl, 1mM DTT, 1mM EDTA, 2mM β -ME, 20% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human NOL3/ARC/NOP (His tag) recombinant protein