

RTN4RL1 Polyclonal Antibody



Catalog Number: 160223

Product name

RTN4RL1 Polyclonal Antibody

Specificity

Mouse, Rat

Antibody description

Polyclonal antibody to RTN4RL1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 220-420 of human RTN4RL1 (NP_848663.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:1000 - 1:4000

Validations



Western blot - RTN4RL1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using RTN4RL1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 20s.