

Prothymosin alpha antibody



Catalog Number: 114223

Product name

Prothymosin alpha antibody

Specificity

Human; other species not tested.

Antibody description

Prothymosin alpha Rabbit Polyclonal antibody.
Positive IHC detected in human lymphoma tissue.

Preparation

This antibody was obtained by immunization of Prothymosin alpha recombinant protein (Accession Number: NM_002823). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions

Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human lymphoma using Catalog No:114223(PTMA antibody) at dilution of 1:100 (under 25x lens)