

Fetuin-A antibody

Catalog Number: 110623

Product name

Fetuin-A antibody

Specificity

Human; other species not tested.

Antibody description

Fetuin-A Rabbit Polyclonal antibody. Positive WB detected in human liver tissue, human heart tissue. Positive IHC detected in human liver tissue. Observed molecular weight by Western-blot: 58 kDa

Preparation

This antibody was obtained by immunization of Fetuin-A recombinant protein (Accession Number: NM_001622). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IHC, WB


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


Validations



human liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:110623(AHSG antibody) at dilution of 1:500



Immunohistochemical of paraffin-embedded human liver using Catalog No:110623(AHSG antibody) at dilution of 1:200 (under 10x lens)



Immunohistochemical of paraffin-embedded human liver using Catalog No:110623(AHSG antibody) at dilution of 1:200 (under 40x lens)