

RFC2 antibody



Catalog Number: 114623

Product name

RFC2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RFC2 Rabbit Polyclonal antibody. Positive WB detected in PC-3 cells, HeLa cells, MCF-7 cells, Raji cells. Observed molecular weight by Western-blot: 35-40 kDa

Preparation

This antibody was obtained by immunization of RFC2 recombinant protein (Accession Number: NM_002914). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:200-1:1000

Validations



PC-3 cells were subjected to SDS PAGE followed by western blot with Catalog No:114623(RFC2 Antibody) at dilution of 1:300