

uPAR antibody

Catalog Number: 116623

Product name

uPAR antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

uPAR Rabbit Polyclonal antibody. Positive IP detected in RAW 264.7 cells. Positive WB detected in RAW 264.7 cells, A2780 cells. Observed molecular weight by Western-blot: 35-70kd

Preparation

This antibody was obtained by immunization of uPAR recombinant protein (Accession Number: XM_047438925). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



RAW264.7 cells were subjected to SDS PAGE followed by western blot with Catalog No:116623(uPAR, PLAUR antibody) at dilution of



1:500

IP Result of anti-uPAR, PLAUR (IP:Catalog No:116623, 4ug; Detection:Catalog No:116623 1:300) with RAW 264.7 cells lysate 4000ug.