

Phospho-MAP3K5-T918 pAb



Catalog Number: 166323

Product name

Phospho-MAP3K5-T918 pAb

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to Phospho-MAP3K5-T918

Preparation

Antigen: A synthetic phosphorylated peptide around T918 of human MAP3K5 .

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB 1:500 - 1:2000

Validations



Western blot - Phospho-MAP3K5-T918 pAb

Western blot analysis of extracts of HeLa cells, using Phospho-MAP3K5-T918 antibody at 1:1000 dilution. HeLa cells were treated by UV for 15-30 minutes. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% BSA.