

Anti-PP-T (2F2F7) Mouse antibody



Catalog Number: 168323

Product name

Anti-PP-T (2F2F7) Mouse antibody

Ig Type

Mouse IgG2a

Specificity

Human

Applications

WB

Antibody description

PP-T (2F2F7) Mouse monoclonal antibody

Dilutions

WB: 1/1000

Preparation

Antigen: Purified recombinant human PP5 protein fragments expressed in E.coli.

Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.


Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Validations



Western blot detection of PP5 in Hela, Jurkat, 293T, A431 and MCF7 cell lysates using PP5 mouse mAb (1:1000 diluted). Predicted band size: 57kDa. Observed band size: 57kDa.