

# Anti-PP-T (2F2F7) Mouse antibody



Catalog Number: 168323

## Product name

Anti-PP-T (2F2F7) Mouse antibody

## Specificity

Human

## Antibody description

PP-T (2F2F7) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human PP5 protein fragments expressed in E.coli.

## Formulation

Ascites containing with 0.03% Proclin300 and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2a


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of PP5 in HeLa, Jurkat, 293T, A431 and MCF7 cell lysates using PP5 mouse mAb (1:1000 diluted). Predicted band size: 57kDa. Observed band size: 57kDa.