

ATPBD4 antibody

Catalog Number: 108323

Product name

ATPBD4 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ATPBD4 Rabbit Polyclonal antibody. Positive IF detected in NIH/3T3 cells. Positive WB detected in NIH/3T3 cells, L02 cells. Observed molecular weight by Western-blot: 48 kDa

Preparation

This antibody was obtained by immunization of ATPBD4 recombinant protein (Accession Number: BC008485). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IF


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

Validations



NIH/3T3 cells were subjected to SDS PAGE followed by western blot with Catalog No:108323(ATPBD4 antibody) at dilution of 1:1000



Immunofluorescent analysis of NIH/3T3 cells using Catalog No:108323(ATPBD4 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)