

VAX2 antibody

Catalog Number: 116723

Product name

VAX2 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

VAX2 Rabbit Polyclonal antibody. Positive IP detected in Y79 cells. Positive WB detected in mouse eye tissue. Observed molecular weight by Western-blot: 31-36 kDa

Preparation

This antibody was obtained by immunization of VAX2 recombinant protein (Accession Number: NM_012476). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IP: 1:200-1:1000

Validations



mouse eye tissue were subjected to SDS PAGE followed by western blot with Catalog No:116723(VAX2 Antibody) at dilution of 1:300



IP Result of anti-VAX2 (IP:Catalog No:116723, 4ug; Detection:Catalog No:116723 1:300) with Y79 cells lysate 2000ug.