

Anti-Glyt2 Rabbit antibody



Catalog Number: 168723

Product name

Anti-Glyt2 Rabbit antibody

Specificity

Human, Mouse, Rat

Antibody description

Glyt2 Rabbit monoclonal antibody

Preparation

Antigen: A synthetic peptide of human Glyt2

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG


Applications

WB

Dilutions

WB: 1/1000

Validations



Western blot detection of Glyt2 in Jurkat,C6 cell lysates using Glyt2 Rabbit mAb(1:1000 diluted).Predicted band size:87kDa.Observed band size:87kDa.