

# DTNA antibody

Catalog Number: 110033

## Product name

DTNA antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

DTNA Rabbit Polyclonal antibody. Positive WB detected in mouse heart tissue, mouse brain tissue, mouse skeletal muscle tissue. Positive IP detected in mouse brain tissue. Observed molecular weight by Western-blot: 49-59 kDa, 78-85 kDa

## Preparation

This antibody was obtained by immunization of DTNA recombinant protein (Accession Number: XM\_047437339). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:1000-1:10000

IP: 1:500-1:5000

## Validations



mouse heart tissue were subjected to SDS PAGE followed by western blot with Catalog No:110033(DTNA Antibody) at dilution of 1:2000



IP Result of anti-DTNA (IP:Catalog No:110033, 4ug; Detection:Catalog No:110033 1:1000) with mouse brain tissue lysate 4000ug.