

NCAPH2 antibody



Catalog Number: 113033

Product name

NCAPH2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

NCAPH2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HEK-293T cells, HepG2 cells. Observed molecular weight by Western-blot: 50 kDa, 70 kDa

Preparation

This antibody was obtained by immunization of NCAPH2 recombinant protein (Accession Number: BC001298). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113033(NCAPH2 Antibody) at dilution of 1:1500