

# ZNF774 antibody

Catalog Number: 117033

## Product name

ZNF774 antibody

## Specificity

Human; other species not tested.

## Antibody description

ZNF774 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 55-63 kDa

## Preparation

This antibody was obtained by immunization of ZNF774 recombinant protein (Accession Number: NM\_001004309). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IP


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

## Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:117033(ZNF774 Antibody) at dilution of 1:600



IP Result of anti-ZNF774 (IP:Catalog No:117033, 3ug; Detection:Catalog No:117033 1:600) with HepG2 cells lysate 2400ug.