

# Anti-Bcl10 Rabbit antibody



Catalog Number: 169033

## Product name

Anti-Bcl10 Rabbit antibody

## Specificity

Human

## Antibody description

Bcl10 Rabbit monoclonal antibody

## Preparation

Antigen: A synthetic peptide of human Bcl10

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Rabbit IgG

## Applications


WB, FC

## Dilutions

1:1000-1:5000

1:50-1:200

## Validations



Western blot detection of Bcl10 in K562,C6,3T3,Hela cell lysates using Bcl10 Rabbit mAb(1:1000 diluted).Predicted band size:26kDa.Observed band size:26kDa.