

F9 antibody



Catalog Number: 110433

Product name

Polyclonal

F9 antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

F9 Rabbit Polyclonal antibody. Positive WB detected in Tunicamycin treated HepG2 cells. Observed molecular weight by Western-blot: 56 kDa

Dilutions


Recommended Dilution:

WB: 1:200-1:2000

Preparation

This antibody was obtained by immunization of F9 recombinant protein (Accession Number: NM_000133). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Tunicamycin treated HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110433(F9 Antibody) at dilution of 1:600