

# F9 antibody

Catalog Number: 110433

---

## Product name

F9 antibody

## Specificity

Human; other species not tested.

## Antibody description

F9 Rabbit Polyclonal antibody. Positive WB detected in Tunicamycin treated HepG2 cells. Observed molecular weight by Western-blot: 56 kDa

## Preparation

This antibody was obtained by immunization of F9 recombinant protein (Accession Number: NM\_000133). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:200-1:2000

## Validations



Tunicamycin treated HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:110433(F9 Antibody) at dilution of 1:600