

TSHR antibody



Catalog Number: 116433

Product name

TSHR antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

TSHR Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, mouse ovary tissue, mouse thymus tissue. Observed molecular weight by Western-blot: 62-70kd

Preparation

This antibody was obtained by immunization of TSHR recombinant protein (Accession Number: NM_001018036). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:116433(TSHR antibody) at dilution of 1:1000