

BCMO1 antibody

Catalog Number: 108433

Product name

BCMO1 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

BCMO1 Rabbit Polyclonal antibody. Positive IHC detected in human testis tissue. Positive WB detected in COLO 320 cells, mouse colon tissue. Observed molecular weight by Western-blot: 65 kDa

Preparation

This antibody was obtained by immunization of BCMO1 recombinant protein (Accession Number: NM_017429). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200

Validations



COLO 320 cells were subjected to SDS PAGE followed by western blot with Catalog No:117106(BCMO1 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:117106(BCMO1 Antibody) at dilution of 1:50 (under 40x lens)