

Human gp130/CD130 (His & Fc Tag) recombinant protein



Catalog Number: 501833

General Information

Gene Name Synonym

CDw130; Interleukin-6 signal transducer; Membrane glycoprotein 130; Oncostatin-M receptor subunit alpha

Protein Construction

A DNA sequence encoding the extracellular domain (Met 1-Ile 618) of human IL6ST (NP_002175.2) precursor was fused with the C-terminal polyhistidine-tagged Fc region of human IgG1 at the C-terminus.

Organism

Human

Expression Host

Human Cells

QC Testing

Activity

1. Measured by its ability to bind human IL11Ra in a functional ELISA.
2. Measured by its ability to inhibit the IL6 R α enhancement of IL6 activity on M1 mouse myeloid leukemia cells (Saito, T. et al. 1991, J. Immunol. 147:168.). The ED₅₀ for this effect is typically 0.05-0.15 μ g/ml in the presence of 50 ng/ml recombinant human IL6sR and 100 ng/ml recombinant human IL6.

Purity

> 95 % as determined by SDS-PAGE

Endotoxin

< 1.0 EU per μ g of the protein as determined by the LAL method

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Glu 23

Molecular Mass

The recombinant human IL6ST/Fc is a disulfide-linked homodimer after removal of the signal peptide. The reduced monomer consists of 844 amino acids and has a predicted molecular mass of 96 kDa. In SDS-PAGE under reducing conditions, the apparent molecular mass of rh IL6ST/Fc monomer is approximately 125-140 kDa due to glycosylation.

Formulation

- Lyophilized from sterile PBS, pH 7.4
1. 5 % trehalose and mannitol are added as protectants before lyophilization.
 2. Please contact us for any concerns or special requirements.

Usage Guide

Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE

Human gp130/CD130 (His & Fc Tag) recombinant protein



Catalog Number: 501833

Human IL6ST / gp130 / CD130 Protein (His & Fc Tag) SDS-PAGE

