

PIR antibody

Catalog Number: 113833

Product name

PIR antibody

Immunogen

[Human PIR Recombinant protein \(GST tag\)](#)

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

PIR Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HepG2 cells, HeLa cells. Observed molecular weight by Western-blot: 35 kDa

Preparation

This antibody was obtained by immunization of PIR recombinant protein (Accession Number: NM_003662). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:1000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:113833(PIR antibody) at dilution of 1:500



IP Result of anti-PIR (IP:Catalog No:113833, 4ug; Detection:Catalog No:113833 1:300) with HeLa cells lysate 1200ug.