

RPS20 antibody



Catalog Number: 114833

Product name

RPS20 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

RPS20 Rabbit Polyclonal antibody. Positive WB detected in MCF7 cells, HEK-293 cells, HepG2 cells, Jurkat cells. Observed molecular weight by Western-blot: 16 kDa

Preparation

This antibody was obtained by immunization of RPS20 recombinant protein (Accession Number: NM_001023). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



MCF7 cells were subjected to SDS PAGE followed by western blot with Catalog No:114833(RPS20 antibody) at dilution of 1:500